

SFWRENG 3K04 – Software Development

Assignment 2

Pacemaker Simulink/DCM

Gurleen Kaur Rahi (400377038)

Yash Bhatia (400362372)

Daniel Lupas (400396848)

Varunan Varathan (400375287)

Connor Usaty (400409624)

Table of Contents

1	Introduction.....	Error! Bookmark not defined.
2	Variables.....	Error! Bookmark not defined.
2.1	Measured Variables	Error! Bookmark not defined.
2.2	Parameters Variables	3
2.3	Controlled Variables.....	6
2.4	Internal Variables.....	7
3	Modules	9
3.1	Pacemaker Modes.....	9
3.1.1	VOO/AOO/VOOR/AOOR	9
3.1.2	VVI/AAI/VVIR/AAIR.....	11
3.1.3	SELECTING PARAMETERS.....	12
3.2	Rate Adaptivity	13
3.2.1	Description.....	13
3.2.2	Variables.....	13
3.2.3	Requirements	14
3.2.4	Anticipated Changes	15
3.3	Hardware Interface	15
3.3.1	PACING	15
3.3.2	SENSING	18
3.4	DCM Communication	19
3.4.1	In.....	19
3.4.2	Out.....	21
4	Design Decisions.....	24
4.1	Pacemaker Modes.....	24
4.1.1	VOO/AOO/VOOR/AOOR	25
4.1.2	AAI/VVI/AAIR/VVIR.....	25
4.2	Rate Adaptivity	26
4.3	Hardware Interface	25
4.3.1	Pacing	29
4.3.2	Sensing	29
4.4	DCM Communication	30
4.4.1	In.....	32
4.4.2	Out.....	33
5	Description of Changes.....	33
6	Testing	34
6.1	VOO.....	34
6.2	AOO.....	35
6.3	VVI.....	35
6.4	AAI.....	37
6.5	AOOR.....	39
6.6	VOOR.....	41
6.7	AAIR.....	42
6.8	VVIR.....	45

1. Introduction

An overview of the Simulink implementation of this Pacemaker Project, written with the end-user in mind, is found in this documentation. This documentation contains the working of ten sensing modes with the goal of creating a functional pacemaker that responds to its surroundings and operates in accordance with those responses.

Depolarization of the ventricles or atria is what is meant by pacing. Sensing is the process of identifying chamber signals. Heart sensing and pacing are provided by the VVI, AAI, VVIR, and AAIR modes. The ventricles employ VVI and VVIR, whereas the atrium uses AAI and AAIR. These four modes operate on an inhibited basis, which means that when a specific activity occurs from the appropriate chamber, the pacemaker is turned off. Only the atrium and/or ventricle are pace (not sensed) by AOO, VOO, AOOR, and VOOR.

While the AOO, VOO, AAI, and VVI modes pace regardless of the patient's activity, the AOOR, VOOR, AAIR, and VVIR modes pace adaptively based on data from the onboard accelerometer.

2. Variables

2.1 Measured Variables

Name	Units/Type	Description	Range
t	ms	Duration of last heartbeat	—
ATR_SIGNAL	Double	Atrial signal assessment on an ECG	0 – 1
VENT_SIGNAL	Double	Ventricle signal assessment on an ECG	0 – 1
ATR_CMP_DETECT	Boolean	Higher than threshold atrial signal voltage	{True, False}
VENT_CMP_DETECT	Boolean	Higher than threshold ventricular signal voltage	{True, False}
Accel	Double x 3	Appropriate pacemaker Accel in g in local (a,b,c)	-4 – 4
Receive	UInt8 x 17	Buffer Accel UART	0 – 255
Buffer_Accel_status	UInt8	shows the fullness of the UART Buffer_Accel	{0, 32}

2.2 Parameters Variables

Name	Units/Type	Description	Range
------	------------	-------------	-------

p_mode	int	Mode of operation in which pacemaker is working	{VOO,AOO,VVI,AAI, DOO, VOOR, AOOR, VVIR, AAIR, DOOR}
p_rateAdaptiveLowrateInterval	ppm	Lowest value of the heart rate that is allowed	30 – 175 ± 8 ms
p_atrialAmplitude	V	Pulse amplitude delivered to atrium	0 – 5 ± 12%
p_ventricularAmplitude	V	Pulse amplitude delivered to ventricle	0 – 5 ± 12%
p_atrial_sensitivity	V	Minimum voltage value that classifies an atrium signal as pulse	0 – 5 ± 2%
p_ventricle_sensitivity	V	Minimum voltage value that classifies a ventricle signal as pulse	0 – 5 ± 2%
p_atrialPulseWidth	ms	Pulse width of atrium	1 – 30 ± 0.2 ms
p_ventricularPulseWidth	ms	Pulse width of ventricle	1 – 30 ± 0.2 ms
p_maximumSensorRate	ppm	Maximum rate of detection indicated by rate adaptivity	50 – 175 ± 4 ppm
p_activityThreshold	–	Rate adaptivity's minimal	{V-Low, Low, Med-Low, Med, Med-High, High, V-high}

		activity level	
p_lowrateLimit	s	Pacing rate's rising time from the lower rate limit to the maximum sensor rate	$10 - 50 \pm 3$ sec
p_response	–	shows reaction to activity levels over the cutoff	{True, False}
p_upperrateLimit	min	Maximum value of the pacing rate decrease that is allowed	$2 - 16 \pm 30$ sec
p_responseFactor	ms	Gives value of general refractory period rather than specific to its type	$150-500 \pm 8$ ms
p_hysterisisLimit	-	bpm interval between the pacemaker's pulse and the last pulse that the body generates on its own	-

2.3 Controlled Variables

Name	Units/Type	Description	Range
PACE_CHARGE_CTRL	Boolean	PWM attached to the primary capacitor	{True, False}

ATR_PACE_CTRL	Boolean	Atrium ring attached to the primary capacitor	{True, False}
VENT_PACE_CTRL	Boolean	Ventricle ring attached to the primary capacitor	{True, False}
ATR_GND_CTRL	Boolean	Atrium ring attached to the ground	{True, False}
VENT_GND_CTRL	Boolean	Ventricle ring connected to ground	{True, False}
Z_ATR_CTRL	Boolean	Circuit with an impedance coupled to the atrium ring	{True, False}
Z_VENT_CTRL	Boolean	Circuit with an impedance coupled to the ventricle ring	{True, False}
PACING_REF_PWM	Percent	Percentage reference PWM needed for primary capacitor	0–100
ATR_CMP_REF_PWM	Percent	For the atrium signal comparator, use reference PWM.	0–100
VENT_CMP_REF_PWM	Percent	For the ventricle signal comparator, use reference PWM.	0–100
s_paceStart	Boolean	Output will be paceStart regardless of Atrium or Ventricle	{True, False}
s_atrialPaceStart	Boolean	Commencement of atrial pulse	{True, False}
s_ventricularPaceStart	Boolean	Commencement of ventricular pulse	{True, False}
s_AOORState	Boolean	Commencement of AOOR state	{True, False}
s_AOORState	Boolean	Commencement of AOO state	{True, False}
s_VOORState	Boolean	Commencement of VOOR state	{True, False}
s_VOORState	Boolean	Commencement of VOO state	{True, False}
s_AAIRState	Boolean	Commencement of AAIR state	{True, False}
s_AAIRState	Boolean	Commencement of AAI state	{True, False}

s_VVIRState	Boolean	Commencement of VVIR state	{True, False}
s_VVIState	Boolean	Commencement of VVI state	{True, False}
s_currentState	Boolean	Refers to the current state	0-5
s_offState	Boolean	States are turned off	{True, False}

2.4 Internal Variables

Name	Units/Type	Description	Range
h_atrial_pulse_detect	Boolean	Pulse detection that happens in atrium	{True, False}
h_ventricle_pulse_detect	Boolean	Pulse detection that happens in ventricle	{True, False}
s_atrialPaceStart	Boolean	Atrium signal pacing gets signalled	{True, False}
s_ventricular_PaceStart	Boolean	Ventricle signal pacing gets signalled	{True, False}
MODES_SENSOR_RATE_Range	ppm	Shows rate adaptivity in the range from minimum to maximum	p_rateAdaptiveLowrateInterval – p_maximumSensorRate

3. Modules

a. Pacemaker Modes

i. VOO/AOO/VOOR/AOOR

1. Description

One chamber is paced at the specified rate in the VOO, AOO, VOOR, and AOOR modes. The different types of variables that were used for these four modes are listed in the table below.

2. Types of Variables

As seen below, we only used measured and internal variables. There were no Programmable Parameters and Controlled Variables used for these modes.

3.1.1.2.1 Measured Variables

Name	Units/Type	Description	Range
t	ms	The duration of time period measured from the previous pulse.	—

3.1.1.2.2 Programmable Parameters

There were no Programmable Parameters used for these modes.

3.1.1.2.3 Controlled Variables

Name	Units/Type	Description	Range
s_AOORState	Boolean	Commencement of AOOR state	{True, False}
s_AOORState	Boolean	Commencement of AOO state	{True, False}
s_VOORState	Boolean	Commencement of VOOR state	{True, False}
s_VOORState	Boolean	Commencement of VOO state	{True, False}

3.1.1.2.4 Internal Variables

Name	Units/Type	Description	Range
AOO_VOO_pace	Boolean	Pacing for either AOO or VOO starts to happen	{True, False}
MODES_SENSOR_RATE	ppm	Indicates the sensor rate for different modes	p_rateAdaptiveLowrateInterval – p_maximumSensorRate
PACE_RATE_INC	ppm/sec	Increase in the pacing rate	0 – 14.5

PACE_RATE_DEC	ppm/sec	Decrease in the pacing rate	0 – 1.3
CURR_PACE_RATE	ppm	Current pacing rate	30 – 175 ppm

3.1.1.3 Requirements

Source	Case	Condition		Activity	End Activity	Destination State	Start Activity
START	N/A	N/A		N/A	N/A	PACING_RATE	CURR_PACE_RATE = MODES_SENSOR_RATE
PACING_RATE	N/A	MODES_SENS OR_RATE > CURR_PACE_R ATE	CURR_PACE_R ATE + PACE_RATE_I NC * t > MODES_SENS OR_RATE	CURR_PACE_RA TE = MODES_SENSOR _RATE	N/A	RATE_UPDATE D	AOO_VOO_pace = True
			CURR_PACE_R ATE + PACE_RATE_I NC * t ≤ MODES_SENS OR_RATE	CURR_PACE_RA TE += PACE_RATE_INC * t			
		MODES_SENS OR_RATE ≤ CURR_PACE_R ATE	CURR_PACE_R ATE – PACE_RATE_D EC * t < MODES_SENS OR_RATE	CURR_PACE_RA TE = MODES_SENSOR _RATE			
			CURR_PACE_R ATE – PACE_RATE_D EC * t ≥ MODES_SENS OR_RATE	CURR_PACE_RA TE -= PACE_RATE_DE C * t			
RATE_UPDATED	N/A	after(60/ CURR_PACE_RATE, sec)		N/A	t = elapsed(sec)	PACING_RATE	on after(1, msec): AOO_VOO_pace = False

3.1.1.4 Anticipated Changes

A possible change that we can do is to combine the logic of these four modes with VVI/AAI/VVIR/AAIR modes to modify the behavior of detecting pulses.

3.1.2 VVI/AAI/VVIR/AAIR

3.1.2.2 Description

In these four modes, pacing occurs in a chamber at a specified rate, however, they are inhibited by the sensed pulses in the same chamber.

3.1.2.3 Types of Variables

As seen below, we only used measured and internal variables. There were no Programmable Parameters and Controlled Variables used for these modes.

3.1.2.3.1 Measured Variables

Name	Units/Type	Description	Range
t	ms	The duration of time period measured from the previous pulse.	—

3.1.2.3.2 Programmable Parameters

There were no Programmable Parameters used for these modes.

3.1.2.3.3 *Controlled Variables*

Name	Units/Type	Description	Range
s_AAIRState	Boolean	Commencement of AAIR state	{True, False}
s_AAIState	Boolean	Commencement of AAI state	{True, False}
s_VVIRState	Boolean	Commencement of VVIR state	{True, False}
s_VVIState	Boolean	Commencement of VVI state	{True, False}

3.1.2.3.4 *Internal Variables*

Name	Units/Type	Description	Range
AAI_VVI_pace	Boolean	The pacing of either AAI or VVI start	{True, False}
MODES_SENSOR_RATE	ppm	Indicates the sensor rate for different modes	p_rateAdaptiveLowrateInterval – p_maximumSensorRate
PACE_RATE_INC	ppm/sec	Increase in the pacing rate	0 – 14.5
PACE_RATE_DEC	ppm/sec	Decrease in the pacing rate	0 – 1.3
Buffer_Accel	double	Last 25 values of Accel	-4 – 4
AV_Pulse_Detect	Boolean	Detecting the pulse generally and not for the type	{True, False}
P_RESPONSEFACTOR	ms	Gives value of general refractory period rather than specific to its type	150–500 ± 8 ms
CURR_PACE_RATE	ppm	Current pacing rate	30 – 175 ppm

3.1.2.4 *Requirements*

Source	Condition			Condition Activities	Finished Activities	Destination State	Starting Activities
START	N/A			N/A	N/A	HOLD	CURR_PACE_RATE = MODES_SENSOR_RATE
HOLD	AV_Pulse_Detect == 1			N/A	t = elapsed(sec)	SENSING_UPDATE	
	after(60/CURR_PACE_RATE – 0.001P_RESPONSEFACTOR + 0.001,sec)			N/A	t = elapsed(sec)	PACING_RATE	
PACING_RATE	after(p_responseFactor, msec)	MODES_SENSOR_RATE >	CURR_PACE_RATE + PACE_RATE_INC	CURR_PACE_RATE =	t += elapsed(sec)	HOLD	AAI_VVI_pace = True on after(1, msec)

		CURR_PACE_RATE	*t > MODES_SENSOR_RATE	MODES_SENSOR_RATE			AAI_VVI_pace = True
			CURR_PACE_RATE + PACE_RATE_INC * t ≤ MODES_SENSOR_RATE	CURR_PACE_RATE + PACE_RATE_INC * t			
		MODES_SENSOR_RATE ≤ CURR_PACE_RATE	CURR_PACE_RATE - PACE_RATE_DEC * t < MODES_SENSOR_RATE	CURR_PACE_RATE = MODES_SENSOR_RATE			
			CURR_PACE_RATE - PACE_RATE_DEC * t ≥ MODES_SENSOR_RATE	CURR_PACE_RATE = PACE_RATE_DEC * t			
SENSING_UPDATE							N/A

3.1.2.5 Possible Changes

For now, there are no possible changes that we can think of doing for these modes.

3.1.3 SELECTING PARAMETERS

3.1.3.1 Description

To feed parameter information to the real mode logic, parameter selection parses it. The mode logic output is then routed to the appropriate signal outputs.

3.1.3.2 Types of Variables

3.1.3.2.1 Measured Variables

No measured variables were used when selecting parameters.

3.1.3.2.2 Programmable Parameters

Name	Units/Type	Description	Range
p_mode	int	Mode of operation in which pacemaker is working	{VOO, AOO, VVI, AAI, DOO, VOOR, AOOR, VVIR, AAIR, DOOR}
p_rateAdaptiveLowrateInterval	ppm	Lowest value of the heart rate that is allowed	30 – 175 ± 8 ms

3.1.3.2.3 Controlled Variables

No Controlled Variables were used when selecting parameters.

3.1.3.2.4 Internal Variables

Name	Units/Type	Description	Range
h_atrial_pulse_detect	Boolean	Pulse detection that happens in atrium	{True, False}
h_ventricle_pulse_detect	Boolean	Pulse detection	{True, False}

		that happens in ventricle	
AOO_VOO_pace	Boolean	The pacing of either AOO or VOO start	{True, False}
AAI_VVI_pace	Boolean	The pacing of either AAI or VVI start	{True, False}
s_atrialPaceStart	Boolean	Atrium signal pacing gets signalled	{True, False}
s_ventricular_PaceStart	Boolean	Ventricle signal pacing gets signalled	{True, False}
AV_CHAMBER	Boolean	Indicates about the type of chamber that's placed next	{A, V}
MODES_SENSOR_RATE_Range	ppm	Shows rate adaptivity in the range from minimum to maximum	p_rateAdaptiveLowrateInterval – p_maximumSensorRate
PACE_RATE_INC	ppm/sec	Increase in the pacing rate	0 – 14.5
PACE_RATE_DEC	ppm/sec	Decrease in the pacing rate	0 – 1.3
AV_Pulse_Detect	Boolean	Detecting the pulse generally and not for the type	{True, False}
P_RESPONSEFACTOR	ms	Gives value of general	150–500 ± 8 ms

		refractory period rather than specific to its type	
--	--	--	--

3.1.3.3 Requirements

Mo des	MODES_SENSOR_RATE	PACE_RATE_INC	PACE_RATE_DEC	AV_Pulse_Detect	P_RESPONSE_FACTOR	S_atrialPaceStart	S_ventricular_PaceStart	AV_CHAMBER
AOO	p_rateAdaptiveLowrateInterval	0	0	X	X	AOO_VOO_pace	0	Atr
VOO				X	X	0	AOO_VOO_pace	Vent
AAI				h_atrial_pulse_detect_Pulse_Detect	p_ap_responseFactor	AAI_VVI_pace	0	Atr
VVI				h_ventricle_pulse_detect_Pulse_Detect	p_vp_responseFactor	0	AAI_VVI_pace	Vent
AOR	MODES_SENSOR_RATE	PACE_RATE_INC	PACE_RATE_DEC	X	X	AOO_VOO_pace	0	Atr
VOO				X	X	0	AOO_VOO_pace	Vent
AAIR				h_atrial_pulse_detect_Pulse_Detect	p_ap_responseFactor	AAI_VVI_pace	0	Atr
VVIR				h_ventricle_pulse_detect_Pulse_Detect	p_vp_responseFactor	0	AAI_VVI_pace	Vent

3.1.3.4 Anticipated Changes

Increase in number of modes would eventually increase the number of programmable parameters so as to accommodate for input and output signals.

3.2 Rate Adaptivity

3.2.1 Description

Based on the activity detected by the on-board accelerometer, the rate adaptivity module provides the appropriate pace. In accordance with the rate adaptivity factors, it also determines the pace at which changes.

3.2.2 Types of Variables

3.2.2.1 Measured Variables

Name	Units/Type	Description	Range
Accel	Double x 3	Appropriate pacemaker Accel in g in local (a,b,c)	-4 – 4

3.2.2.2 Programmable Parameters

Name	Units/Type	Description	Range
p_rateAdaptiveLowrateInterval	ppm	Lowest value of the heart rate that is allowed	30 – 175 ± 8 ms
p_maximumSensorRate	ppm	Maximum rate of detection indicated by rate adaptivity	50 – 175 ± 4 ppm

p_activityThreshold	–	Rate adaptivity's minimal activity level	{V-Low, Low, Med-Low, Med, Med-High, High, V-high}
p_lowrateLimit	s	Pacing rate's rising time from the lower rate limit to the maximum sensor rate	10 – 50 ± 3 sec
p_response	–	shows reaction to activity levels over the cutoff	{True, False}
p_upperrateLimit	min	Maximum value of the pacing rate decrease that is allowed	2 – 16 ± 30 sec

3.2.2.3 Controlled Variables

There are no Controlled Variables for Rate Adaptivity.

3.2.2.4 Internal Variables

Name	Units/Type	Description	Range
PACE_RATE_INC	ppm/sec	Increase in the pacing rate	0 – 14.5
PACE_RATE_DEC	ppm/sec	Decrease in the pacing rate	0 – 1.3
MODES_SENSOR_RATE	ppm	Indicates the sensor rate for different modes	p_rateAdaptiveLowrateInterval – p_maximumSensorRate
Buffer_Accel	double	Last 25 values of Accel	-4 – 4
PEAK_TO_PEAK_SIG	double	Peak to peak of activity signal	0 – 4

3.2.3 Requirements

Variable	Value
PEAK_TO_PEAK_SIG	max(Buffer_Accel) – min(Buffer_Accel)
PACE_RATE_INC	$(p_maximumSensorRate - p_rateAdaptiveLowrateInterval) / p_lowrateLimit$
PACE_RATE_DEC	$(p_maximumSensorRate - p_rateAdaptiveLowrateInterval) / (60 * p_upperrateLimit)$
Overdriving_Action	$4 * p_response * (p_PEAK_TO_PEAK_SIG - p_activityThreshold)$
MODES_SENSOR_RATE	$\min(p_maximumSensorRate, \max(p_rateAdaptiveLowrateInterval, Overdriving_Action))$

3.2.4 Possible Changes

Change in the design for rate adaptivity requires modifying this document a lot as all the previous requirements will be changed.

3.3 Hardware Interface

3.3.1 Pacing

3.3.1.1 Description

In addition to controlling which level to charge the pacing capacitor at and which chamber to pace, the pacing module also keeps an eye on the pace start and AV choose lines.

3.3.1.2 Types of Variables

3.3.1.2.1 Measured Variables

There were no measured variables used in this section.

3.3.1.2.2 Programmable Parameters

Name	Units/Type	Description	Range
p_atrialAmplitude	V	Pulse amplitude delivered to atrium	0 – 5 ± 12%
p_ventricularAmplitude	V	Pulse amplitude delivered to ventricle	0 – 5 ± 12%
p_atrialPulseWidth	ms	Pulse width of atrium	1 – 30 ± 0.2 ms
p_ventricularPulseWidth	ms	Pulse width of ventricle	1 – 30 ± 0.2 ms

3.3.1.2.3 Controlled Variables

Name	Units/Type	Description	Range
PACE_CHARGE_CTRL	Boolean	PWM attached to the primary capacitor	{True, False}
ATR_PACE_CTRL	Boolean	Atrium ring attached to the primary capacitor	{True, False}
VENT_PACE_CTRL	Boolean	Ventricle ring attached to the primary capacitor	{True, False}
ATR_GND_CTRL	Boolean	Atrium ring attached to the ground	{True, False}
VENT_GND_CTRL	Boolean	Ventricle ring connected to ground	{True, False}
Z_ATR_CTRL	Boolean	Circuit with an impedance coupled to the atrium ring	{True, False}
Z_VENT_CTRL	Boolean	Circuit with an impedance coupled to the ventricle ring	{True, False}
PACING_REF_PWM	Percent	Percentage reference PWM needed for primary capacitor	0–100

3.3.1.2.4 Internal Variables

Name	Units/Type	Description	Range
s_atrialPaceStart	Boolean	Commencement of atrial pulse	{True, False}
s_ventricular_PaceStart	Boolean	Commencement of ventricular pulse	{True, False}
AV_CHAMBER	Boolean	Indicates about the type of chamber that's placed next	{A, V}

3.3.1.3 Requirements

Source	Case	Condition	Condition Activities	Finished Activities	Destination State	Starting Activities
START	N/A	N/A	N/A	N/A	Y_c22Charged	VENT_GND_CTRL=False Z_VENT_CTRL=False Z_ATR_CTRL=False ATR_GND_CTRL=False
Y_c22Charged	N/A	s_ventricular_PaceStart == False AND s_atrialPaceStart == False	N/A	N/A	Buffer_Accel	VENT_PACE_CTRL=False ATR_PACE_CTRL=False PACE_CHARGE_CTRL=True AV_BLOCK_CTRL=True
Buffer_Accel	N/A	s_atrialPaceStart == True	N/A	N/A	Y_atriumPaced	N/A
		s_ventricular_PaceStart == True	N/A	N/A	Y_ventriclePaced	N/A
y_atriumPaced	N/A	after(p_atrialPulseWidth)	N/A	N/A	Discharge_21_Atr	PACE_CHARGE_CTRL=False VENT_PACE_CTRL=False VENT_GND_CTRL=False ATR_GND_CTRL=False ATR_PACE_CTRL=True
Y_initialDischarge	N/A	N/A	N/A	N/A	Y_c22Charged	ATR_PACE_CTRL=False ATR_GND_CTRL=True
Y_ventriclePaced	N/A	after(p_ventricularPulseWidth)	N/A	N/A	Discharge_21	PACE_CHARGE_CTRL=False ATR_PACE_CTRL=False ATR_GND_CTRL=False VENT_GND_CTRL=False VENT_PACE_CTRL=True
Y_initialDischarge	N/A	N/A	N/A	N/A	Y_c22Charged	VENT_PACE_CTRL=False; VENT_GND_CTRL=True;

3.3.1.4 Possible Changes

Some modifications in configuration will occur when any change in hardware will be made.

3.3.2 Sensing

3.3.2.1 Description

When natural cardiac pulses are recognized, the sensor module analyzes the signal, sets the values of the reference capacitors, and outputs pulses.

3.3.2.2 Types of Variables

3.3.2.2.1 Measured Variables

Name	Units/Type	Description	Range
ATR_CMP_DETECT	Boolean	Higher than threshold atrial signal voltage	{True, False}
VENT_CMP_DETECT	Boolean	Higher than threshold ventricular signal voltage	{True, False}

3.3.2.2.2 Programmable Parameters

Name	Units/Type	Description	Range
p_atrial_sensitivity	V	Minimum voltage value that classifies an atrium signal as pulse	0 – 5 ± 2%
p_ventricle_sensitivity	V	Minimum voltage value that classifies a ventricle signal as pulse	0 – 5 ± 2%

3.3.2.2.3 Controlled Variables

Name	Units/Type	Description	Range
ATR_CMP_REF_PWM	Percent	For the atrium signal comparator, use reference PWM.	0–100
VENT_CMP_REF_PWM	Percent	For the ventricle signal comparator, use reference PWM.	0–100
LEAD_SENSE_CIRCUIT	Boolean	Lead-connected sensing circuit	{True, False}

3.3.2.2.4 Internal Variables

Name	Units/Type	Description	Range
h_atrial_pulse_detect	Boolean	Pulse detection that happens in atrium	{True, False}
h_ventricle_pulse_detect	Boolean	Pulse detection that happens in ventricle	{True, False}

3.3.2.3 Requirements

Variable	Value
h_atrial_pulse_detect	ATR_CMP_DETECT
h_ventricle_pulse_detect	VENT_CMP_DETECT
ATR_CMP_REF_PWM	p_atrial_sensitivity / (5 V) * 100 %
VENT_CMP_REF_PWM	p_ventricle_sensitivity / (5 V) * 100 %
LEAD_SENSE_CIRCUIT	1

3.3.2.4 Possible Changes

Configuration will have to be modified if the any change is made in the hardware design.

3.4 DCM Communication

3.4.1 Data Input

3.4.1.1 Description

In order to update the pertinent variables on the pacemaker, the communication-Data Input module watches the UART Receive signal for messages from the DCM.

3.4.1.2 *Types of Variables*

3.4.1.2.1 *Measured Variables*

Name	Units/Type	Description	Range
Receive	Uint8 x 17	Buffer_Accel UART	0 – 255
Buffer_Accel_status	Uint8	shows the fullness of the UART Buffer_Accel	{0, 32}

3.4.1.2.2 *Programmable Parameters*

Name	Units/Type	Description	Range
p_mode	int	Mode of operation in which pacemaker is working	{VOO, AOO, VVI, AAI, DOO, VOOR, AOOR, VVIR, AAIR, DOOR}
p_rateAdaptiveLowrateInterval	ppm	Lowest value of the heart rate that is allowed	30 – 175 ± 8 ms
p_AV_delay	ms	Delay that is caused between ventricle and atrial signals	70 – 300 ms ± 8 ms
p_atrialAmplitude	V	Pulse amplitude delivered to atrium	0 – 5 ± 12%
p_ventricularAmplitude	V	Pulse amplitude delivered to ventricle	0 – 5 ± 12%
p_atrial_sensitivity	V	Minimum voltage value that classifies an atrium signal as pulse	0 – 5 ± 2%
p_ventricle_sensitivity	V	Minimum voltage value that classifies a ventricle signal as pulse	0 – 5 ± 2%
p_atrialPulseWidth	ms	Pulse width of atrium	1 – 30 ± 0.2 ms
p_ventricularPulseWidth	ms	Pulse width of ventricle	1 – 30 ± 0.2 ms
p_vp_responseFactor	ms	Refractory time period of ventricle	150–500 ± 8 ms
p_ap_responseFactor	ms	Refractory time period of atrium	150–500 ± 8 ms
p_maximumSensorRate	ppm	Maximum rate of detection indicated by rate adaptivity	50 – 175 ± 4 ppm

p_activityThreshold	–	Rate adaptivity's minimal activity level	{V-Low, Low, Med-Low, Med, Med-High, High, V-high}
p_lowrateLimit	s	Pacing rate's rising time from the lower rate limit to the maximum sensor rate	10 – 50 ± 3 sec
p_response	–	shows reaction to activity levels over the cutoff	{True, False}
p_upperrateLimit	min	Maximum value of the pacing rate decrease that is allowed	2 – 16 ± 30 sec

3.4.1.2.3 Controlled Variables

There were no controlled variables used in this section.

3.4.1.2.4 Internal Variables

Name	Units/Type	Description	Range
Received_Code	uint8	Last received Received_Code from DCM	1– 5

3.4.1.3 Requirements

Buffer_Accel_status	Receive(1)	Parameters	Received_Code
0	1	Receive(2:17)	Receive(1)
	X	no change	
32	X	no change	X

3.4.1.4 Possible Changes

Reconfiguration will be required if any change in the hardware is made as the requirements would be different.

3.4.2 Data Output

3.4.2.1 Description

To transmit the pacemaker's serial number, echo parameter data, or atrial/ventricular signal data to the DCM, the communications-Data output module regulates the UART UART_BUFFER_ACCEL signal.

3.4.2.2 Variables

3.4.2.2.1 Measured Variables

Name	Units/Type	Description	Range
ATR_SIGNAL	Double	Atrial signal assessment on an ECG	0 – 1

VENT_SIGNAL	Double	Ventricle signal assessment on an ECG	0 – 1
-------------	--------	---------------------------------------	-------

3.4.2.2.2 Programmable Parameters

Name	Units/Type	Description	Range
p_mode	int	Mode of operation in which pacemaker is working	{VOO,AOO,VVI,AAI, DOO, VOOR, AOOR, VVIR, AAIR, DOOR}
p_rateAdaptiveLowrateInterval	ppm	Lowest value of the heart rate that is allowed	30 – 175 ± 8 ms
p_AV_delay	ms	Delay that is caused between ventricle and atrial signals	70 – 300 ms ± 8 ms
p_atrialAmplitude	V	Pulse amplitude delivered to atrium	0 – 5 ± 12%
p_ventricularAmplitude	V	Pulse amplitude delivered to ventricle	0 – 5 ± 12%
p_atrial_sensitivity	V	Minimum voltage value that classifies an atrium signal as pulse	0 – 5 ± 2%
p_ventricle_sensitivity	V	Minimum voltage value that classifies a ventricle signal as pulse	0 – 5 ± 2%
p_atrialPulseWidth	ms	Pulse width of atrium	1 – 30 ± 0.2 ms
p_ventricularPulseWidth	ms	Pulse width of ventricle	1 – 30 ± 0.2 ms
p_vp_responseFactor	ms	Refractory time period of ventricle	150–500 ± 8 ms
p_ap_responseFactor	ms	Refractory time period of atrium	150–500 ± 8 ms
p_maximumSensorRate	ppm	Maximum rate of detection indicated by rate adaptivity	50 – 175 ± 4 ppm
p_activityThreshold	–	Rate adaptivity's minimal activity level	{V-Low, Low, Med-Low, Med, Med-High, High, V-high}
p_lowrateLimit	sec	Pacing rate's rising time from the lower rate limit to the	10 – 50 ± 3 sec

		maximum sensor rate	
p_response	–	shows reaction to activity levels over the cutoff	{True, False}
p_upperrateLimit	min	Maximum value of the pacing rate decrease that is allowed	2 – 16 ± 30 sec

3.4.2.2.3 Controlled Variables

Name	Units/Type	Description	Range
UART_BUFFER_ACCEL	Uint8 x 16	UART escape Buffer_Accel	0– 255

3.4.2.2.4 Internal Variables

Name	Units/Type	Description	Range
Received_Code	uint8	Last received Received_Code from DCM	1– 5

3.4.2.3 Requirements

Received Code	UART_BUFFER_ACCEL
0	Serial Digit
1	X
2	Parameters
3	[ATR_SIGNAL, 0]
4	[VENT_SIGNAL, 0]
5	[ATR_SIGNAL, VENT_SIGNAL]

3.4.2.4 Possible Changes

Reconfiguration will be required if any change in the hardware is made as the requirements would be different.

4 Design Decisions

4.1 Pacemaker Modes




Figure 1: Pacing Modes Subsystem

The pacemaker's primary logic is contained in three Stateflow charts, which mostly interact with abstracted variables to successfully accomplish hardware concealing. For AOO, VOO, AOOR, and VOOR modes, there is a single Stateflow chart; for AAI, VVI, AAIR, and VVIR modes, but only one is active at any given moment. With the exception of input and output signal routing (which is managed by 4.1.4 Parameter Selection), all of the grouped modes are implemented in essentially the same way. Because Received_Code is reduced, sharing Stateflow charts facilitates improved maintainability.

In every mode, when pacing, the current pacing rate is adjusted at the appropriate rate based on the duration of the previous waiting interval if the rate reported by the sensor differs from the current pacing rate. It is assumed that the rate of change is linear in time. The higher rate limit was determined to have

no purpose in these modes and was therefore omitted from this design as maximum sensor rate is defined as being independently programmable from the upper rate limit. On the appropriate out signal, a 1 ms pulse indicates that pacing should start. The hardware interface is informed of the expected level of the next pacing by the AV choose signal.

4.1.1 AOO/VOO/AOOR/VOOR




Figure 2: AOO/VOO/AOOR/VOOR Stateflow chart

Since they don't need to sense the heart, the AOO, VOO, AOOR, and VOOR modes just alternate between pacing and non-pacing states.

4.1.2 AAI/VVI/AAIR/VVIR




Figure 3: AAI/VVI/AAIR/VVIR Stateflow chart

The XOO modes' structure is comparable to that of the AAI, VVI, AAIR, and VVIR modes. Additionally, if a spontaneous pulse is detected after the refractory period (during which spontaneous pulses are ineffective), it is possible to avoid the pacing condition. To make sure the pacemaker does not pace

4.2 Rate Adaptivity



24

It is expected that the pace of increase and decline in pace rate is linear in time. For comprehensive equations, see 4.2.3.

4.3 Hardware Interface




Figure 5a: Hardware Interface subsystem




Figure 6b: Hardware Interface Input Variables




Figure 7c: Hardware Interface Output Pins

The pacemaker's IN and OUT variables are abstracted by the hardware interface as a response factor, making it simple to transfer the pacemaker's primary functionality to any hardware device. This guarantees that the logic describing how the chosen instructions should be carried out on the hardware would be the only modifications that are required. In order to output 1 ms pulses when atrial/ventricular signals surpass their respective sensitivities, the hardware interface only has to accept 1 ms pulsed pace atrium/ventricle signals and an atrial-ventricular level selector signal.

4.3.1 Pacing




Figure 8: Pacing Stateflow chart

When a pulse is detected on the S_atrialPaceStart or S_ventricular_PaceStart lines, the pacing module detects it and waits in a Buffer_Accel state. Then, it sets the relevant flags to pace the heart. Flags are set to start charging the pacing capacitor again and discharging the blocking capacitor after a delay equal to the pulse width. As a precaution against looping caused by a non-deasserted pacing signal pulse, the Buffer_Accel state is used.

4.3.2 Sensing




Figure 9: Sensing Subsystem

To identify natural pulses, the sensing module routes the flags to the variables h_atrial_pulse_detect and h_ventricular_pulse_detect, and sets the proper PWMs for the atrial and ventricular reference capacitors.

4.4 DCM Communication




Figure 10a: Communications Subsystems




Figure 11b: Communications Subsystems for DCM Communication

For all external communication, the DCM Communication module is responsible. Like the hardware interface, it interacts with the particular IN and OUT variables while presenting abstracted variables to the system. This design selected UART communication.

4.4.1 In




Figure 12: Communications In Stateflow chart

When the Buffer_Accel_status indicates that it is full, the In module reads the Receive Buffer_Accel and, according on the Received_Code received, either updates the parameter data or transmits data to the DCM.

4.4.2 Out




Figure 13: Communications Out Stateflow chart

Only upon receiving a Received_Code of a number between 1 and 16 (both inclusive) does the out module become active. As described in 4.4.2.3, data is transmitted via UART.

5 Description of Changes

A. We performed the following actions to add Rate adaptivity:

1. With the same inputs and outputs, every previously implemented mode (AOO, VOO, AAI, and VVI) was retained. The new input source for the lowrateLimit is the only difference.
2. Another subsystem that was in charge of implementing Rate Adaptivity received the lowrateLimit. The subsystem did not alter the lowrateLimit whether the modes were AXX or VXX. The Rate Adaptivity algorithm was used to adjust the lowrateLimit based on Accelerometer data if the modes were AXXR or VXXR.
3. We used a linear function to apply the acceleration to a pulse in bpm. For this algorithm, the maximum value was the maximumSensorRate, and the minimum was the lowrateLimit.
4. ActivityThreshold, Reaction Time, Response Factor, and Recovery Time were included as input variables.

B. We did the following to introduce hysteresis:

1. This modification had no effect on AOO/AOOR or VOO/VOOR.
2. The AAI/AAIR and VVI/VVIR subsystems saw a total modification in the stateflow.
3. The system detects if hysteresis was turned on when it detects a "sensed" or natural pulse. If so, it would not wait for the lowrateLimit time, but for the hysteresisLimit time. If the pacemaker had only pulsated, this would not occur because in that case the wait time would always be lowrateLimit.
4. Every output is same. It was necessary to add more inputs to the AAI/AAIR and VVI/VVIR modes. Hysteresis and its Limitation.

C. We carried out the following actions to add UART Communications:

1. We removed every single constant block that controlled our variables. It was possible to create a UART system that output variables as UINT8 or UINT16 variables. This subsystem was in charge of sending egram data to the dcm, making adjustments to the pacemaker, and repeating data from the pacemaker back to the dcm.
2. Before being reinserted into the variables of the current system, the outputs were transformed to the appropriate units and data types.
3. This subsystem needed the following inputs: the device id, VENT_SIGNAL, and ATR_SIGNAL. Each and every output had customizable parameters.

6 Testing


6.1 AOO

Chamber paced: Atrium
Chamber sensed: N/A
Response to Sensing: N/A

Test Cases:

- a) Natural Atrium: OFF | Natural Ventricle: OFF

Only one test case—in which the ventricle and atrium are both nonfunctional—can show that our pacemaker is functioning correctly in AOO mode because we are not sensing any chambers.



Graph representing artificial pulse from HeartView

6.2 VOO

Chamber paced: Ventricle


Chamber sensed: N/A

Response to Sensing: N/A

Test Cases:

a) Natural Atrium: OFF | Natural Ventricle: OFF

Only one test case—in which the ventricle and atrium are both nonfunctional—can show that our pacemaker is functioning correctly in VOO mode as well because we are not sensing any chambers.



Graph representing artificial pulse from HeartView

6.3 AAI

Chamber paced: Atrium

Chamber sensed: Atrium


Response to Sensing: Inhibited

Test Cases:

a) Natural Atrium: ON, Pulse Width: 1ms | Natural Ventricle: OFF | Heart Rate: 30bpm

Because the heart beats every 2000 milliseconds and the pacemaker is designed to maintain a heart rate of 60 beats per minute, which is a pulse after every 1000


milliseconds, the pacemaker is expected to create a pulse at very low heart rate and low pulse width. The graph below illustrates how our pacemaker fills the gap by sending a pulse every 1000 ms.



Graph representing natural and artificial pulse from HeartView

- b) Natural Atrium: ON, Pulse Width: 10ms | Natural Ventricle: OFF | Heart Rate: 30bpm


The pacemaker is expected to generate a pulse when the pulse width is increased to 10 ms while maintaining the same heart rate. This is because there is a gap of more than 1000 ms between two consecutive pulses. The graph below illustrates how our pacemaker accomplishes the same thing by supplying the pulse needed to fill the gap in time.



Graph representing natural and artificial pulse from HeartView

- c) Natural Atrium: ON, Pulse Width: 1ms | Natural Ventricle: OFF | Heart Rate: 50bpm


A pacemaker is intended to produce an artificial pulse when our heart rate is slightly below the natural rate but the pulse width is insufficient to achieve the natural rate, which is 60 beats per minute or one pulse every 1000 milliseconds. It is evident from the graph below that our pacemaker compensates for the delay.



Graph representing natural and artificial pulse from HeartView

- d) Natural Atrium: ON, Pulse Width: 10ms | Natural Ventricle: OFF | Heart Rate: 60bpm

A pacemaker should not produce an artificial pulse when the heart rate is below normal but the pulse width is sufficient to fill the gap. Our pacemaker isn't giving our heart an extra pulse, as the graph below illustrates.



Graph representing natural pulse from HeartView

6.4 VVI


Chamber paced: Ventricle
Chamber sensed: Ventricle
Response to Sensing: Inhibited

Test Cases:

- a) Natural Atrium: OFF | Natural Ventricle: ON, Pulse Width: 1ms | Heart Rate: 30bpm

The pacemaker's pulse should be generated frequently enough to sustain a normal heart rate even at very low heart rates.


The graph below shows that our pacemaker generates artificial pulses when it has been more than 1000 milliseconds since the last pulse, which aids in sustaining a normal heart rate of 60 beats per minute or a pulse every 1000 milliseconds.



Graph representing natural and artificial pulse from HeartView

- b) Natural Atrium: OFF | Natural Ventricle: ON, Pulse Width: 10ms | Heart Rate: 30bpm


High pulse width cannot fill the void left by a low heart rate; therefore, a pacemaker is needed to ensure that the heart is beating properly. It is evident from the graph below that our pacemaker keeps the heart operating normally by giving it artificial pulses.



Graph representing natural and artificial pulse from HeartView

- c) Natural Atrium: OFF | Natural Ventricle: ON, Pulse Width: 1ms | Heart Rate: 50bpm

A pacemaker is supposed to activate and start sending pulses to the heart when the heart is beating somewhat slower than usual but the pulse width is not wide enough to make up for the difference. We may observe from the graph below that our pacemaker performs as planned.




Graph representing natural and artificial pulse from HeartView

- d) Natural Atrium: **OFF** | Natural Ventricle: **ON**, Pulse Width: **10ms** | Heart Rate: **60bpm**

A pacemaker shouldn't be used when the patient's heart rate is slightly below normal and their pulse width is sufficient to close the gap.

The graph below shows that our pacemaker is not causing the heart to beat artificially.



Graph representing natural pulse from HeartView

6.5 AOOR


Chamber paced: Atrium

Chamber sensed: N/A


Response to Sensing: N/A

Rate Modulation

Heartview signals: Natural Atrium: **OFF** | Natural Ventricle: **OFF**




a) No Activity sensed



Pace to Atrium with an LRL of 60 ppm if no activity is detected above the threshold.

b) Activity sensed



Increase LRL within a 10-second response time and a 2-minute recuperation time after detecting Activity above Threshold.

Result: Passed

6.6 VOOR

Chamber paced: Ventricle

Chamber sensed: N/A


Response to Sensing: N/A

Rate Modulation

Heartview signals: Natural Atrium: **OFF** | Natural Ventricle: **OFF**




a) No Activity sensed



Pace to Ventricle with an LRL of 60 ppm if no activity is detected above the threshold.

b) Activity Sensed



Increase LRL after 10 seconds of detecting Activity above Threshold and within 2 minutes of recovering from it.

Result: Passed

6.7 AAIR

Chamber paced: Atrium

Chamber sensed: Atrium

Response to Sensing: Inhibited


Rate Modulation

Test Cases:


a) No activity

Heartview signals: Natural Atrium: **On** | PW: 1ms | Heart rate: 30Bpm

DCM signals:



The pacemaker's job is to fill the gap in heart rate while we are in AAIR mode, with no activity, and when the heart rate is lower than normal. Our pacemaker functions in the same way, as seen by the graph below.



Graph representing artificial and natural pulse from HeartView

b) Physical activity

Heartview signals: Natural Atrium: **On** | PW: 1ms | Heart rate: 30Bpm

DCM signals:

Welcome to the Device Control Module

CONNECTED

Disconnect
Change Mode
Show Egram

Current Device UID:
0xF00D


Logged In As:
admin

Current Device Mode:
AAIR

Current Parameters

Lower Rate Limit: 60	Edit LRL
Upper Rate Limit: 120	Edit URL
Atrial Amplitude: 5.0	Edit AAMP
Atrial PW: 1	Edit APW
ARP: 250	Edit ARP
Atrial Sensitivity: 5.0	Edit ASens
PVARP: 250	Edit PVARP
Hysteresis: OFF	Toggle
Max Sensor Rate: 120	Edit Max Rate
Activity Threshold: 0	Update AT
Reaction Time: 30	Update RT
Response Factor: 8	Update RF
Recovery Time: 5	Update RT

Compared to normal, our pacemaker beats more frequently while we are physically active.




Graph representing artificial and natural pulse from HeartView


c) No activity

Heartview signals: Natural Atrium: **On** | PW: 1ms | Heart rate: 60Bpm

DCM signals:



A pacemaker is not meant to provide an extra pulse while the heart is beating naturally and there is no physical activity. The graph below makes it evident that our pacemaker functions in the same way.




Graph representing artificial and natural pulse from HeartView


d) **Physical activity**

Heartview signals: Natural Atrium: **On** | PW: 1ms | Heart rate: 60Bpm

DCM signals:



The purpose of a pacemaker is to add a second pulse to account for physical activity when the patient's heart is beating regularly. The graph below makes it evident that our pacemaker functions the same way.



Graph representing artificial and natural pulse from HeartView

6.8 VVIR


Chamber paced: Ventricle
 Chamber sensed: Ventricle
 Response to Sensing: Inhibited
 Rate Modulation

Test Cases:


a) No activity

Heartview signals: Natural Ventricle: **On** | PW: 1ms | Heart rate: 30Bpm

DCM signals:



The pacemaker's job is to provide a pulse to fill the gap when it is in the VVIR mode, there is no activity, and the heart rate is lower than normal. The graph below shows that our pacemaker functions in the same way.




Graph representing artificial and natural pulse from HeartView


b) Physical activity

Heartview signals: Natural Ventricle: **On** | PW: 1ms | Heart rate: 30Bpm

DCM signals:



Our pacemaker beats more frequently during physical activity than it does during rest.




Graph representing artificial and natural pulse from HeartView


c) No activity

Heartview signals: Natural Ventricle: **On** | PW: 1ms | Heart rate: 60Bpm

DCM signals:



A pacemaker is not meant to provide an extra pulse while the heart is beating naturally and there is no physical activity. The graph below makes it evident that our pacemaker functions in the same way.



Graph representing artificial and natural pulse from HeartView